

NOW READY
TO USE

500 IOT
2000
3000 ml
PERFADEX® PLUS

1500
1000
Perfadex® Plus is a pre-buffered, extracellular solution of dextran 40 containing calcium. The device is intended for the flushing, cold static storage and transportation of isolated lungs after removal from the donor in preparation for eventual transplantation into a recipient.

FR : Perfadex® Plus est une solution extra-cellulaire pré-tamponnée à base de Dextran 40, contenant du calcium. Le dispositif est indiqué pour le rinçage, le stockage statique à froid et le transport de poumons isolés après retrait sur le donneur et en préparation à une éventuelle transplantation chez un receveur.

REF 99450



DEX

STERILE



25 ± 25°C



500 Rx Only
Manufactured by
XVIVO Perfusion AB
Milestone data 10.

PERFADEX® Plus

THE GOLD STANDARD IN LUNG
PRESERVATION MADE READY TO USE

Combines safety and efficacy of PERFADEX® with easier use

XVIVO
PERFUSION

PERFADEX® Plus

THE GOLD STANDARD IN LUNG PRESERVATION MADE READY TO USE

PERFADEX® Plus – a new solution for optimal cold static preservation of donor lungs

XVIVO Perfusion have developed PERFADEX® Plus a new generation of PERFADEX®. To save you time and effort, we have pre-supplemented the solution with calcium ions and THAM, so it is always ready to use. Calcium ions have been added to further mimic the composition of plasma.



05

New easy to use click port – with one simple click, the XVIVO Click Adapter connects to the bag. Once connected, the tubing can rotate on its own axis for flexibility during handling.

PERFADEX® Plus has several design features that make it even safer and easier to handle.

01

Strong, reinforced hanger – for safety and convenience.

02

Free from PVC, latex and phthalates - for patient safety.

03

New double, twist-off sterile ports are easy to grip and easy to open.

04

Injection port is made of stronger plastic – no risk of perforating the bag during injection.

PERFADEX® Plus is an extracellular, low potassium, dextran-based electrolyte preservation solution for rapid cooling, perfusion and cold static storage of donor lungs **pre supplemented with calcium ions and THAM so that it is ready to use whenever you need it.**

How PERFADEX® Plus works

PERFADEX® Plus is indicated for the flushing, cold static storage and transportation of isolated lungs after removal from the donor in preparation for eventual transplantation into a recipient. The colloid component, dextran 40, protects the microvasculature against post-ischemic reperfusion injury, primarily by preventing pathological leukocyte-endothelial interaction^{12,13} It also prevents edema formation during preservation.^{14,15,20}

Numerous studies have shown that PERFADEX® enables safe preservation of lungs up to 12 hours, depending on the quality of the donor lung being retrieved.^{1,2,3,4,5,6,7,8,9,10,11}

Using PERFADEX® Plus with click port

PERFADEX® Plus must be stored between 2-25 °C. The solution must not be supplemented with calcium ions or THAM – PERFADEX® Plus is pre-supplemented so that it is ready to use whenever you need it.

PERFADEX® Plus now comes in a bag that allows for easier handling - with only one click the bag is connected to XVIVO Click Adapter which gives you full flexibility during the perfusion process as it can rotate around its own axis.

Once opened the solution should be kept chilled and used within 24 hours.

Current best practice²¹

According to Munshi et al, 2013

Current best practice based on the latest research on potential lung donor and lung-preservation techniques advises the use of an extracellular, dextran-based solution.

Lung preservation techniques

- ✓ Extracellular solution consisting of dextran-40, glucose and low potassium
- ✓ Anterograde and retrograde flushing of 60 mL/kg and max 30 cm height
- ✓ Temperature during lung preservation 2-8 °C
- ✓ Inflation to 50% of total lung capacity, fraction of inspired oxygen 50%
- ✓ Pharmacological additives: prostaglandin E1, heparin, glucocorticoids
- ✓ Cold ischemic times generally less than 8 h
- ✓ Normothermic ex-vivo lung perfusion based on lung assessment and therapeutics.

The Endothelium - a vulnerable tissue

The lung is primarily composed of endothelial cells which line the enormous surface area of the capillaries (equivalent to an entire tennis court) and a similar surface area of types I and II epithelial cells which line the alveoli and secrete surfactant respectively. The endothelium is the most vulnerable and plays a critical role for the structure and function of a normal vessel wall. Endothelial cells produce a variety of biologically active substances that control vascular permeability, vessel tone, coagulation, fibrinolysis and inflammatory responses. Some of these substances, such as proteins which seal the junctions between cells (adhesion molecules), are integral parts of the cell structure. Others, such as nitric oxide (NO), prostacycline, chemokines, or factors involved in coagulation and fibrinolysis, are produced and then released by the endothelial cells either lumenally or abuminally.¹⁶

Consequences of an injured endothelium

Injuries to the endothelium can induce platelet and leukocyte sticking. This can trigger inflammatory cascades including increased permeability of the capillary wall, which in turn increases tissue edema and the risk of Primary Graft Dysfunction (PGD). A well preserved endothelium is antithrombogenic, yet promotes platelet aggregation and coagulation if injured.^{17,18}


The importance of an intact endothelium

Experimental and clinical evidence indicates that early ischemia-reperfusion injury to the endothelium, within the very first few hours of reperfusion, is a key trigger for cytokine cascades which eventually lead to PGD and subsequent graft failure often months or years later.¹⁸ This early injury can be prevented or mitigated by minimizing physical injury in manipulation and storing the lungs in a protective solution under optimal temperature conditions.¹⁹

What causes damage to the endothelium?

A number of factors can injure the pulmonary endothelium during the manipulation and temporary storage involved in the retrieval of donor lungs;

- Traumatic manipulation during retrieval, evaluation and transplantation
- Excessive pressure
- Low temperature – particularly below 2°C
- Storage solution – e.g. intracellular type (high K+) solutions
- Prolonged cold ischemia
- Ischemia-reperfusion – free radical injury



Excess oncotic pressure retains water in capillaries and prevents loss into the extravascular (interstitial) space, thus preventing lung edema. The coating of the endothelium with dextran 40 prevents excessive leukocyte interaction.

US - approved by the FDA

Usage: PERFADEX® Plus must be stored between 2-25°C (36-77°F). Must not be supplemented with calcium ions or THAM – PERFADEX® Plus is pre-supplemented so that it is ready to use whenever you need it. PERFADEX® Plus is supplied in 1000ml and 3000ml bags.

Intended use: PERFADEX® Plus solution for lung perfusion is indicated for flushing, cold static storage and transportation of isolated lungs after removal from the donor in preparation for eventual transplantation into a recipient.

References

1. Muller C, et al, transplantation, 1999, 68(8):1139-43.
2. Fischer S, et al, J Thorac. Cardiovasc Surg, 2001, 121(3):594-6
3. Struber M, et al, Eur. J. Cardiothorac. Surg, 2001, 19 (2):190-4.
4. Rabanal J M, et al, Transplant Proc, 2003, 35(5):1938-9
5. Sakamaki F, et al, Am J Respir Crit Care Med, 1997;156:1073-81
6. Struber M, et al, J Heart Lung Transplant, 1999;18:187
7. Struber M, et al, J Thorac Cardiovasc Surg, 2000, 120 (3):566-72.
8. Oto T et al, Ann Thorac Surg, 2006 Nov;82(5):1842-8.
9. Gamez P et al Arch Bronconeumol. 2005 Jan;41(1):16-19.
10. Okada Y et al, Ann Thorac Cardiovasc Surg, 2006 Feb;12(1):10-4
11. Amaoutakis GJ et al, J Heart Lung Transplant, 2010 Dec; 29(12): 1380-7
12. Menger M D, Transplant Proceed, 1995, 27; 5:2863-65
13. Hoffman h et al, Z. Hertz-, Thorax-, Gefässchir, 1997, 11: 108-14
14. Frost -Arner L et al, Microsurgery, 1995, 16: 357-361
15. Keshavjee SH et al, J Thorac Cardiovascular Surg, 1992; 103:314-25
16. Toborek M, J Kaisers S, Basic Res Cardiol, 1999; 94: 295-314
17. Zilla P et al, J Card Surg, 1993; 8: 32-60
18. Davis S.F. et al, Circulation, 1996; 93: 457-462
19. Steen S, Scand Cardiovasc J, 2001; 35: 297-301
20. Arfors K-E, Buckley, PB, Bailliere's Clinical Anaesthesiol, 1997; 11, 1, 15-47
21. Munshi L et al, Lancet Respir Med 2013; 1; 318-28

Supply information:

Reference:	Packaging:
19850	10 x 1000 ml bags with PERFADEX® Plus
19950	2 x 3000 ml bags with PERFADEX® Plus