en: Instructions for Use

XIVO Lung Cannula Set™
For connection of isolated donor lungs to an extracorporeal perfusion system.**Description**
The XIVO Lung Cannula Set™ is a sterile, single-use set to be used to connect isolated donor lungs to an extracorporeal perfusion system for ex-vivo assessment.**Intended Use**
The XIVO Lung Cannula Set™ is a sterile, single-use set intended to be used to connect isolated donor lungs to an extracorporeal perfusion system for ex-vivo assessment.

The XIVO Lung Cannula Set™ is not intended for any other use than that stated in the section Intended use.

Clinical Benefit
XIVO Lung Cannula Set™ enables safe Ex Vivo Lung Perfusion (EVLP). EVLP increases the number of donor lungs available for transplant.**Contraindications**
None known.**Storage conditions**
Store between 10–25 °C.**Instructions for Use**
Read and understand all of the instructions below before proceeding. Ensure aseptic technique is undertaken when required.**Straight Cannula (1) Yellow Cap:**
1. Place the straight cannula into the pulmonary artery up to the point of bifurcation.

2. Suture left atrial cuff around cone using 4-0 monofilament.

3. Connect the pressure line to sensor via yellow-capped luer end.

Cone Cannula (1) Green Cap:

1. Cut cones to fit inside left atrial cuff radius.

2. Suture left atrial cuff around cone using 4-0 monofilament.

3. Connect the pressure line to sensor via green-capped luer end.

Concave Cannula (2) White Cap:

1. The concave cannula with the white luer-capped end is to be used as a back-up cannula to either the pulmonary artery or left atrial cuff.

2. Cutt cannula according to fit into the pulmonary artery or left atrial cuff.

3. Tie (or suture) cannula in place.

4. Connect the pressure line to sensor via white-capped luer end.

Dispose the XIVO Lung Cannula Set™ after the procedure according to the hospital's standard routines.

Rx Only - This device is restricted to sale by or on the order of a physician.

Serious incidents information
Note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**Specifications**Sterility (EO) SAL 10⁴

Bakteriell endotoksin [EU/device] <20

Symboler

	Sterilized using Ethylene Oxide		Date of manufacture		Medical Device
	Catalogue number		Single sterile barrier		
	Batch code		Manufacturer		
	Use by date		Importer		
	Store between +10°C and +25°C		Consult Instructions for use		Authorized representative in Switzerland

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