

XVIVO SILICONE TUBING SET™



For perfusion during donor lung procurement

The Silicone Tubing Set™ is specifically designed for rapid flushing of donor lungs during procurement. The set permits simultaneous/continuous perfusion from two PERFADEX® bags or two PERFADEX® Plus bags. Given either together or one at the time, whereby one bag can be exchanged while the other is still flushing.

The wide bore 3/8" tubing permits use of high flow rates at relatively low pressures, in order to minimise the risk for endothelial injury. The use of silicone tubing reduces the risks of kinking.

The tubing set for lung perfusion is 200cm long and consists of two 50cm tubes which are joined together by a y-connector to a 150cm tube. The 50cm tubes are equipped with spikes and all three tubes are fitted with tubing clamps for regulating the flow. The tubing set is sterile and intended for single use only.

XVIVO Silicone Tubing Set™

- Rapid flushing during procurement
- Simultaneous and continuous perfusion from two PERFADEX® bags or two PERFADEX® Plus bags
- High flow rates at relatively low pressures
- Able to regulate flow rates by use of tubing clamps

Supply information:

Reference: 19005	Packaging: 1 tubing set / pack
---------------------	-----------------------------------
